

# AA ATOMIC ABSORPTION SPECTROPHOTOMETRY DATA SHEET No. 9

## Calibration Curves of Standard Solutions Using AA-6200

Shimadzu AA-6200 offers high sensitivity analyses

### Measurement of **Cu**

**Cu**

#### Lamp Condition

Lamp Current : 6mA  
Set Wavelength : 324.8nm  
Slit Width : 0.7nm  
Lamp Mode : BGC-D2


$$\text{Abs.} = 0.17615\text{Conc} + 0.00986338$$

$$r = 0.9998$$

#### Atomizer Condition

Fuel Gas Flow Rate : 1.8 l/min.  
Support Gas Flow Rate : 8.0 l/min.  
Flame Type : Air-C<sub>2</sub>H<sub>2</sub>

Conc.	Abs.
0.0000	0.0004
0.5000	0.0987
1.0000	0.1934
2.0000	0.3670
5.0000	0.8871



### Measurement of **Cr**

**Cr**

#### Lamp Condition

Lamp Current : 10mA  
Set Wavelength : 357.9nm  
Slit Width : 0.7nm  
Lamp Mode : NON-BGC


$$\text{Abs.} = 0.0751206\text{Conc} + 0.000965948$$

$$r = 1.0000$$

#### Atomizer Condition

Fuel Gas Flow Rate : 3.0 l/min.  
Support Gas Flow Rate : 7.2 l/min.  
Flame Type : Air-C<sub>2</sub>H<sub>2</sub>

Conc.	Abs.
0.0000	-0.0002
0.5000	0.0387
1.0000	0.0773
2.0000	0.1512
5.0000	0.3763



### Measurement of **Fe**

**Fe**

#### Lamp Condition

Lamp Current : 12mA  
Set Wavelength : 248.3nm  
Slit Width : 0.2nm  
Lamp Mode : BGC-D2


$$\text{Abs.} = 0.114959\text{Conc} + 0.00376523$$

$$r = 0.9998$$

#### Atomizer Condition

Fuel Gas Flow Rate : 2.2 l/min.  
Support Gas Flow Rate : 8.0 l/min.  
Flame Type : Air-C<sub>2</sub>H<sub>2</sub>

Conc.	Abs.
0.0000	-0.0025
0.5000	0.0621
1.0000	0.1234
2.0000	0.2366
5.0000	0.5764



### Measurement of **Pb**

**Pb**

#### Lamp Condition

Lamp Current : 12mA  
Set Wavelength : 217.0nm  
Slit Width : 0.7nm  
Lamp Mode : BGC-D2


$$\text{Abs.} = 0.0682955\text{Conc} + 0.00385065$$

$$r = 0.9999$$

#### Atomizer Condition

Fuel Gas Flow Rate : 2.0 l/min.  
Support Gas Flow Rate : 7.0 l/min.  
Flame Type : Air-C<sub>2</sub>H<sub>2</sub>

Conc.	Abs.
0.0000	-0.0004
1.0000	0.0735
2.0000	0.1435
5.0000	0.3464
10.0000	0.6855



## AA-6200 Hardware Specifications

<b>Optics</b>	Double Beam (chopper mirror)
<b>Monochromator</b>	aberration corrected Czerny-Turner monochromator Holographic grating
<b>Wavelength range</b>	190-900nm (Automated wavelength selection)
<b>Slit</b>	0.2nm, 0.7nm (Manual setting)
<b>Background correction</b>	D2-Lamp method
<b>Lamp turret</b>	2-lamps (Manual turret)
<b>Nebulizer</b>	Nebulizer integrated impact bead and jacket tip Pt/Ir capillary
<b>Chamber</b>	Polypropylene
<b>Burner</b>	Titanium 10cm slot burner (Optional high temperature burner) Fixed back/forward position and burner height (Simply switching of Air/C <sub>2</sub> H <sub>2</sub> and N <sub>2</sub> O/C <sub>2</sub> H <sub>2</sub> burner)
<b>Gas control</b>	Manual setting of flow rate Automatic Air/N <sub>2</sub> O switching system
<b>Ignition</b>	Push ignite button
<b>Safety</b>	Gas pressure monitoring to prevent flashback Automatic flame monitoring Safety interlock for misuse of burner The flame goes out automatically in case of power failure.
<b>Software</b>	Software based on MS-Windows™ 3.1 or MS-Windows 95™
<b>Dimensions and weight</b>	W690 × D425 × H370mm, 38kg
<b>Power requirements</b>	AC220, 230V ± 10%, 50/60Hz (Certification of CE mark)

## AA-6200 Software Specifications

<b>Data processing</b>	
① Calibration curve	Linear regression analysis by least square method, Abs=f (Conc.) Dimension : 1st, 2nd, 3rd order, ON/OFF of zero intercept Selectable calibration curve by Curve ID
② Accuracy enhancement	Usage of average value, cancel of abnormal data by the limitation of CV or SD values.
③ Sensitivity enhancement	Automatic calculation by sensitivity correction
④ Actual concentration calculation	Automatic calculation of actual concentration by sample volume, fill-up volume, and dilution factor upon data acquisition
<b>Signal process</b>	
① Measurement procedure setting	Easy setting by Wizard function
② Data display	Measurement procedure, analysis results, actual concentration, calculation
③ MRT worksheet	Factor of actual concentration, combined display of data and time display of multiple elements by tab sheet
④ Signal display	Signal profile display, Current signal profile display, Overlay display of signal profiles
<b>Retrieval and storage of data</b>	<b>Automatic loading of template parameter file upon software activation</b> <b>Automatic data storage function</b> <b>Raw signal profile can be stored.</b>
<b>Import/Export of data</b>	Import sample name from ASCII file, Export data file by ASCII file



SHIMADZU CORPORATION. International Marketing Division  
3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan  
Phone: 81(3)3219-5641 Fax. 81(3)3219-5710 Cable Add.:SHIMADZU TOKYO

SHIMADZU SCIENTIFIC INSTRUMENTS, INC.  
7102 Riverwood Drive, Columbia, Maryland 21046, U.S.A.  
Phone: 1(410)381-1227 Fax. 1(410)381-1222 Toll Free: 1(800)477-1227

SHIMADZU DEUTSCHLAND GmbH  
Albert-Hahn-Strasse 6-10, D-47269 Duisburg, F.R. Germany Phone: 49(203)7687-0 Fax. 49(203)766625

SHIMADZU (ASIA PACIFIC) PTE LTD.  
16 Science Park Drive #01-01 Singapore Science Park, Singapore 118227, Republic of Singapore  
Phone: 65-778 6280 Fax. 65-779 2935

SHIMADZU SCIENTIFIC INSTRUMENTS (OCEANIA) PTY. LTD.  
Units F, 10-16 South Street Rydalmere N.S.W. 2116, Australia  
Phone: 61(2)9684-4200 Fax. 61(2)9684-4055

SHIMADZU DO BRASIL COMERCIO LTDA.  
Rua Cenzo Sbrighi, 25, Agua Branca, Sao Paulo, CEP 05036-010, BRAZIL  
Phone: (55)11-3611-1688 Fax. (55)11-3611-2209

SHIMADZU (HONG KONG) LIMITED  
Suite 1028 Ocean Center, Harbour City, Tsim Sha Tsui, Kowloon HONG KONG  
Phone: (852)2375-4979 Fax. (852)2199-7438

**Overseas Offices**  
Istanbul, Beijing, Shanghai, Guangzhou, Shenyang, Chengdu, Moscow

**URL** <http://www.shimadzu.com>